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1. Total Suspended Solids

Definition

Total Suspended Solids are particles that are larger than 2 microns found in the water column, can be trapped by a filter that is analyzed using a filtration apparatus. And then dry the filter and the residue retained on it to constant weight in a 103-105 °C oven.

Principle

It is a water quality parameter used to assess the quality of a specimen of any type of water or water body, ocean water for example, or wastewater after treatment in a wastewater treatment plant. A measured volume of well-mixed sample is vacuum filtered through a dried pre-weighed glass fiber filter. The filters and residue is dried to a constant weight at 103-105°C. The increase in weight of the filter represents the total suspended solids.

2. Volatile Suspended Solids

Definition

Volatile Suspended Solids are a water quality measure obtained from the loss on ignition of the mass of measured total suspended solids. This ignition generally takes place in an oven at a temperature of 550 ± 50 °C in a muffle furnace.

Principle

The residue from the above methods is ignited to constant weight at 550°C. The remaining solids represent the fixed total, dissolved or suspended solids while the weight lost on ignition is the volatile solids. The determination is useful in control of wastewater treatment plant operation because it offers a rough estimate of the amount of organic matter present in the solid fraction of wastewater, activated sludge and industrial wastes.

3. Total Solids (TS)

Definition

The material left in a sample vessel after evaporation and subsequent oven drying at 103-105 °C.

Principle

An appropriate volume of well-mixed sample is evaporated to dryness in a weighed evaporating dish in a hot air oven at 103-105°C, the residue remaining are cooled and weighed. The residual material in the evaporating dish is classified as total solids, and may consist of organic, inorganic, dissolved, suspended or volatile matter.

4. Total Volatile Solids (TVS)

Definition

The total solids lost from a sample after ignition for specified time at 550 °C

Principle

The residue from total solids is ignited to constant weight at 550 ± 50 °C more than 15 minutes in a muffle furnace. The remaining solids are fixed total solids, while those lost to ignition are volatile total solids. The determination is useful in control of wastewater treatment plant operation because it offers a rough estimate of the amount of organic matter present in the solid fraction of wastewater, activated sludge and industrial wastes.

5. Moisture Content

Definition

The moisture content is amount of water consisted in sludge forms which can be evaporated after drying in a 103-105 °C.

Principle

The moisture content is measured by determination of percent of total solid (%TS). The residue from sludge is dried to at constant weight 103-105 °C oven more than 1 hour. And calculate with in % Total solids as below equation;

$$\% \text{ Total Solids} = \frac{(A - B) \times 100}{C - B}$$

Where

A	= final weight of dried residue + dish (mg)
B	= weight of dish (mg)
C	= weight of wet sample + dish (mg)

NATS LAB