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## Definition and Principle of Microbiology Parameters

### 1. Total coliform, Fecal coliform and *E. coli*

#### Definition

Coliform bacteria have been used for many years to determine the quality and safety of water for human consumption. *Escherichia coli* and other groups of coliforms may be present where there has been fecal contamination originating from warm-blooded animals which can be contaminated in various water sources to the environment.

#### Significance

The most probable number (MPN) of bacteria present can then be estimated from the number of tubes inoculated and the number of positive tubes obtained in the confirmatory test, using specially devised statistical tables. It is applicable to the analysis of water of drinking water quality, salt or brackish water as well as muds, sediments, and fecal sludge.

Colony forming unit (CFU) is a method which measures the number of viable bacterial cell in sample. The method which counts colony forming units is referred as standard plate count. The viable colonies that appear on the agar plates are expressed as CFU per 1 ml (colony forming unit per milliliter) of the sample for liquids or CFU per 1 g (colony forming unit per one gram) of the sample for solids.

### 2. Male-specific (F<sup>+</sup>) Coliphage (Bacteriophage MS2)

#### Definition

Coliphage are viruses (bacteriophage) that infect *E. coli* and are indicators of fecal contamination. This method can detect male-specific (F<sup>+</sup>) Coliphage. The adoption of this method standard describes the detection and quantification of MS2 Coliphage in wastewater or effluent, fecal sludge and dried or composted sludge.

### Significance

Quantitation of Bacteriophage, use the specific host for coliphage and was activated by TSA with antibiotics and then incubate at 37 degree Celsius, 24 Hr. After that, Count the plaque that appear on the agar.

## **3. Helminth Egg and Nematode**

### Definition

The viable helminth ova count is indication of the effectiveness of the disinfection process of the wastewater treatment facility and generate data on the prevalence of the different species, as well as provide some indication of the incidence of infection in communities serviced.

The helminth ova and nematode count. NATS laboratory can detection, identification and quantification of total and viable helminth ova in wastewater or effluent, wet sludge or stool sample, faecal sludge and dried or composted sludge.

### Significance

Helminth Egg and Nematode use the different method from characteristic of sample and count by microscope. If it was detected and then incubate in 0.1 N H<sub>2</sub>SO<sub>4</sub> with simulate intestine of humans for verify wastewater treatment system can remove Helminth egg and nematode or not.

## **4. *Clostridium perfringens* (*C. perfringens*)**

### Definition

*Clostridium perfringens* is widely recognized as a valuable indicator for faecal and food pollution. Intestinal tract of animals and man, these Gram-positive bacteria from spores which are resistant to heating compared with vegetative cells. Monitoring of *C. perfringens* has proven useful for the assessment of the quality of water resources include fecal sludge and to check the stages of water treatment to evaluate the treatment-works performance.

### Significance

The quantitation analysis of *C. perfringens*, analyzed by characteristic of sample. Solid sample can observe the gas generation in LS medium and liquid sample can count the black colony on TSC agar.

## 5. DNA Extraction

### Definition

This method can isolate microbial genomic DNA from even tough compost, sediment, manure, biomass and fecal sludge from water treatment system samples for study the biodiversity.

### Significance

DNA extraction is based on the destruction of cell walls using lysozyme or detergent to break cells first. Then digest the proteins and RNA using protease and RNase. Final step is precipitate the DNA using cool absolute ethanol, which separates the DNA into a single strand and will be analyzed by Bioinformatic tool further.

NATS LAB